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Species	Cells	Enzyme(s)	Reference
Human	脂肪 细胞	2型胶原 酶:0.01-0.5%	Effect of Collagenase Concentration on The Isolation of Small Adipocytes from Human Buccal Fat Pad., J Oral Sci, 2018
Human	脂肪 细胞	Collagenase Type 1: 0.1%	Epigenome-wide Association Study of Body Mass Index, and the Adverse Outcomes of Adiposity, Nature541, 81, 2017
Human	基质 细胞	胶原酶2 型:0.075%	Ultrasound-Assisted Liposuction Provides a Source for Functional Adipose-Derived Stromal Cells., Cytotherapy 19, 1491-1500, 2017
Human	骨髓 基质 细胞	胶原酶1 型:0.075%	Mesenchymal Stromal Cells Protect Human Cardiomyocytes from Amyloid Fibril Damage, Cytotherapy 19, 1426-1437, 2017
Human	间充 质干 细胞	胶原酶2 型:0.2%	Adipogenic Differentiation of Mesenchymal Stem Cells Alters Their Immunomodulatory Properties in a Tissue-Specific Manner., Stem Cells 35, 1636-1646, 2017
Human	脂肪 源内 皮细 胞	胶原酶1 型:0.1%	Autologous Cell Sources in Therapeutic Vasculogenesis: In Vitro and In Vivo Comparison of Endothelial Colony-Forming Cells from Peripheral Blood and Endothelial Cells Isolated from Adipose Tissue., Cytotherapy 18, 242-52, 2016
Human	脂肪 干细 胞	胶原酶1 型:0.1%	High Glucose-Induced Reactive Oxygen Species Generation Promotes Stemness in Human Adipose-Derived Stem Cells, Cytotherapy 18, 371-83, 2016
Human	骨髓 基质 细胞	胶原酶4 型:0.2%	Effect of Mild Heat Stress on the Proliferative and Differentiative Ability of Human Mesenchymal Stromal Cells., Cytotherapy17, 359-68, 2015
Human	脂肪 干细 胞	胶原酶2 型:0.1%	Expression Analysis of Human Adipose-Derived Stem Cells During In Vitro Differentiation to an Adipocyte Lineage., BMC Med Genomics 8, 41, 2015
Human	脂肪 基质	胶原酶1 型:0.075%	Therapeutic Potential of Adipose-Derived SSEA-3-Positive Muse Cells for Treating Diabetic Skin Ulcers., Stem Cells Transl Med 4, 146, 2015
Human	脂肪 组织 细胞	胶原酶1 型:0.1%	Differential Effects of Processing Time and Duration of Collagenase Digestion on Human and Murine Fat Grafts., Plast Reconstr Surg 136, 189e-199e, 2015

Human	脂肪 基质 血管 细胞	中性蛋白 酶:2.4 u/ml	Human White and Brite Adipogenesis is Supported by MSCA1 and is impaired by Immune Cells., Stem Cells 33, 1277-91, 2015
Human	脂肪 源性 间充 质干 细胞	胶原酶2 型:0.1%	Defined Serum-Free Media for In Vitro Expansion of Adipose-Derived Mesenchymal Stem Cells., Cytotherapy 16, 915, 2014
Human	脂肪 提取 干细 胞	胶原酶1 型:0.15%	Choosing the Right Type of Serum for Different Applications of Human Adipose Tissue-Derived Stem Cells: Influence on Proliferation and Differentiation Abilities., Cytotherapy 16, 789, 2014
Human	间充 质基 质	胶原酶1 型:0.1%	Proliferative and Phenotypical Characteristics of Human Adipose Tissue-Derived Stem Cells: Comparison of Ficoll Gradient Centrifugation and Red Blood Cell Lysis Buffer Treatment Purification Methods., Cytotherapy 16, 1220-8, 2014
Human	脂肪 间质 干细 胞	动物游离胶 原酶:200u /ml	Xenofree Enzymatic Products for the Isolation of Human Adipose-Derived Stromal/Stem Cells., Tiss Eng 19, 473-8, 2013
Human	基质 血管 成分	胶原酶1 型:0.075%	Stromal Vascular Fraction Isolated from Lipo-Aspirates Using an Automated Processing System: Bench and Bed Analysis., J Tissue Eng Regen Med 7, 864, 2013
Human	脂肪 干细 胞	胶原酶1 型:0.1%	Platelet-Rich Plasma Greatly Potentiates Insulin-Induced Adipogenic Differentiation of Human Adipose-Derived Stem Cells Through a Serine/Threonine Kinase Akt-dependent Mechanism and Promotes Clinical Fat Graft Maintenance., Stem Cells Transl Med 1, 206-20, 2012
Human	血管 周围 干细 胞	胶原酶2 型:0.1%	An Abundant Perivascular Source of Stem Cells for Bone Tissue Engineering., Stem Cells Transl Med 1, 673, 2012
Human	脂肪 来源 的间 质血 管	胶原酶1 型:0.1%	Concise Review: Adipose-Derived Stromal Vascular Fraction Cells and Platelet-Rich Plasma: Basic and Clinical Implications for Tissue Engineering Therapies in Regenerative Surgery., Stem Cells Transl Med 1, 230-6, 2012

Mouse	脂肪 基质	胶原酶1 型:0.1%	Adipose Stromal Vascular Fraction-Mediated Improvements at Late-Stage Disease in a Murine Model of Multiple Sclerosis., Stem Cells 35, 532-544, 2017
Mouse	脂肪 间质 细胞	胶原酶2 型:0.2%	Adipose Mesenchymal Stromal Cells Minimize and Repair Radiation-Induced Oral Mucositis., Cytotherapy 18, 1129-45, 2016
Mouse	脂肪 基质	胶原酶1 型:0.1%	Improved Mobilization of Exogenous Mesenchymal Stem Cells to Bone for Fracture Healing and Sex Difference., Stem Cells 34, 2587-2600, 2016
Mouse	脂肪 细胞	胶原酶:0.1%	The Effects of A Single Developmentally-Entrained Pulse of Testosterone in Female Neonatal Mice On Reproductive and Metabolic Functions in Adult Life., Endocrinology 156, 3737, 2015
Mouse	间质 血管 细胞	胶原酶2 型:0.2%	Natural Killer T Cells in Adipose Tissue are Activated in Lean Mice., Exp Anim 62, 319, 2013
Mouse	脂肪 干细 胞	胶原酶2 型:0.1%	Biological and Clinical Availability of Adipose-Derived Stem Cells for Pelvic Dead Space Repair., Stem Cells Transl Med 1, 803, 2012

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