
GU101019





conjugated)

021-34695924 Order

orders@ab-mart.com

400-6123-828 Support

support1@ab-mart.com

Web www.ab-mart.com.cn

Description:

☐ 5ml;10ml;50ml

Whole IgG antibodies are isolated as intact molecules from antisera by immunoaffinity chromatography. They have an Fc portion and two antigen binding Fab portions joined together by disulfide bonds and therefore they are divalent. The average molecular weight is reported to be about 160 kDa. The whole IgG form of antibodies is suitable for the majority of immunodetection procedures and is the most cost effective. Horseradish peroxidase (HRP) conjugates are prepared by a modified Nakane and Kawaoi procedure (J. Histochem. Cytochem. 1974. 22, 1084). Peroxidase conjugates are commonly used for immunohistochemistry, Western blotting, and ELISA. Affinitypurified anti-horseradish peroxidase and conjugates are available for detection of horseradish peroxidase antigen or for signal amplification of HRP-containing reagents. For immunostaining of mammalian cells, an advantage of using anti-horseradish peroxidase is reduced background, since the antibody does not recognize the endogenous peroxidase-like enzymes found in those cells.

UltraPolymer Goat anti-Rabbit IgG (H&L) conjugated to HRP, affinity purified, min x w/ human, mouse or rat serum.

Species reactivity: Rabbit

Antibody type: Goat polyclonal Antibody

Isotype: Goat IgG Storage buffer:

TBS, 1% (w/v) BSA, 0.1% Proclin 150, 10% heat-inactivated animal serum.

Form: Liquid

Concentration: Ready to Use

Purity: Affinity purified using solid phase Rabbit IgG. Recommended dilutions: IHC-P: Ready to Use

Store at 2-8

Storage: °C.